




DRAINAGE TUBE CODE 110



This product has been designed for clearing bloody, serous or purulent fluid or air collections from body compartments as a result of post-surgical waste (abdominal, pelvic, retroperitoneal surgery).


- Improved patient comfort.
- 450 mm long.
- Versión soft:
 - Enhanced flexibility and reduced stiffness allow catheter to be “milked”.
 - Improved patient comfort.
 - Lengths of 50 and 100 cm without perforations.



Product
100%
Silicone



Open, rounded end.



5 side holes.
Free of latex

full-length radiopaque line.

Code	Size (FR)	Outer Ø (mm)	Inner Ø (mm)
110-12	12,0	4,00	2,50
110-16	16,0	5,30	3,50
110-20	20,0	6,70	4,50
110-24	24,0	8,00	5,50
110-28	28,0	9,30	6,50
110-32	32,0	10,70	7,50
110-36	36,0	12,00	8,50
110-40	40,0	13,30	9,50

Code	Size (FR)	Ø Outer (mm)	Ø Inner (mm)	Length (cm)
110-15	15,0	5,00	3,50	50,00
110-15	15,0	5,00	3,50	100,00
110-18	18,0	6,00	4,00	50,00
110-18	18,0	6,00	4,00	100,00
110-24	24,0	8,00	5,50	50,00
110-24	24,0	8,00	5,50	100,00
110-36	36,0	12,00	8,50	50,00
110-36	36,0	12,00	8,50	100,00

Soft version