




BALLOON GASTROSTOMY TUBE

CODE 330-331



This product has been designed to administer food directly into the gastric cavity in patients unable to feed by mouth, patients with gastroesophageal reflux and upper gastrointestinal tract malformations, gastroesophageal surgery patients and temporary or permanent gastrostomy patients.

- Fixation ring that prevents the displacement of the probe.
- Ball located 1.5 mm from the distal end that remains adhered to the tubulant while it is not inflated.
- The balloon is sealed at the distal end and at the proximal end it has an automatic shut-off valve that allows it to be inflated.
- Connections type Toomey, luer slip and luer lock.



100% silicone with radiopaque line.
Fixation ring prevents tube movement.
abdominal wall.



Wide range of sizes available.



Centimetre marks to measure the abdominal wall.

Paediatric			
Code	Size (FR)	Outer Ø (mm)	Balloon capacity (ml)
331-14	14,0	4,70	5
331-16	16,0	5,30	20

2 ways - 3 accesses

Adult			
Code	Size (FR)	Outer Ø (mm)	Balloon capacity (ml)
330-18	18,0	6,00	20
330-20	20,0	6,70	20
330-22	22,0	7,30	20
330-24	24,0	8,00	20

2 ways - 3 accesses